

Estimation of Alcohol

Advance Chemistry Lab 4

Assist. Lect. Zahra A. Ismail

Alcohols

are organic compounds characterized by the presence of a hydroxyl group in their compounds. They are derived from hydrocarbons by replacing the hydrogen atom with a hydroxyl group. Alcohols have the general formula ROH.



This is the functional group of an alcohol

Principle of Alcohol Estimation

- Alcohol reacts with oxidizing agent i.e. potassium dichromate and gets oxidized to acetic acid.
- Remaining of the potassium dichromate will react with iodine. This is a light sensitive process, so incubation is carried out in dark.
- Liberated iodine is titrated with 0.1N sodium thiosulphate.
- In this process starch is used as an indicator.

Reaction



Experimental Work

Chemicals and equipments

- Potassium Dichromate ($\text{K}_2\text{Cr}_2\text{O}_7$)
- Potassium Iodide (KI) [25%]
- Sulphuric Acid (H_2SO_4)
- Starch [1%]
- Distilled Water
- Sodium Thiosulphate ($\text{Na}_2\text{S}_2\text{O}_3$) [0.1N]
- Flask
- Graduated cylinder
- Pipet


Experiment Steps

- Add 5 ml Distilled water to 1 ml of Alcohol.
- Add 10ml of $\text{K}_2\text{Cr}_2\text{O}_7$ (1 g $\text{K}_2\text{Cr}_2\text{O}_7$ dissolved in 450 ml H_2SO_4) and incubate at room temperature for 10 minute (color is yellow).
- Empty the tube containing Alcohol + $\text{K}_2\text{Cr}_2\text{O}_7$ to fresh flask containing 100ml Distilled water.
- Add 4ml KI (25%) + 2-3 drops of Starch (color becomes blue).
- Add dropwise Sodium thiosulphate (24.8 g Sodium thiosulphate dissolved in 1L water) [0.1N] until the blue color disappears and note the reading.


1gm $\text{K}_2\text{Cr}_2\text{O}_7$ in 450 mL conc. H_2SO_4
Oxidized $\text{K}_2\text{Cr}_2\text{O}_7$



25% Potassium Iodide



1% Starch



24.8 gm in 1 litre water
0.1 N Sodium thiosulfate

