



# Radiation physics –practical

## Experimental No.1

# Linear and mass absorption coefficients

Asst.Lec.Dumoaa Haider Shakir

## 1.1 The objectives of this experiment:

- ▶ 1) Verification of the absorption law of Gamma radiation.
- ▶ 2) Determination:
  - a) The linear absorption coefficient ( $\mu$ ).
  - b) The mass absorption coefficient  $\mu_m$  .
  - c) The half value thickness of the absorbing material ( $X_{1/2}$ )

## 1.2 Apparatus

- ▶ GM Tube and stand (Counter box, power supply – transformer, GM Tube, shelf stand, USB cable, and a source holder for the stand) – shown in Figure 1.
- ▶ Source of radiation.
- ▶ Sheets of different absorbing materials (Aluminum and Lead).



# 1.3 Theory

- ▶ When Gamma radiation passes through matter, it undergoes attenuation primarily by Compton, photoelectric and pair production interactions. The intensity of the radiation is thus decreased as a function of thickness of the absorbing material. The mathematical expression for intensity ( I ) is given by the following expression:

$$I = I_0 e^{-\mu x} \dots \dots \dots (1)$$

where,

$I_0$  the original intensity of the beam.

$I$  is the intensity transmitted through an absorber thickness  $X$ .

$\mu$  is the linear absorption coefficient for the absorbing material

- ▶ If we rearrange eq.(1) and take the natural logarithm of both sides, the expression becomes,

$$\ln \frac{I}{I_0} = -\mu X \dots\dots (2)$$

- ▶ The half value layer (HVL) of the absorbing material is defined as that thickness  $X_{1/2}$  which will decrease the initial intensity by half. That is,  $I = I_0/2$ . If we substitute this into eq.(2) we get:


$$\ln (2) = \mu X_{1/2} \dots\dots (3)$$

- ▶ Rearranging eq.(3) we get:

$$X_{1/2} = \ln (2) / \mu \dots\dots (4) \text{ Or } X_{1/2} = 0.693 / \mu$$

▶ The factors which effects on the Linear and mass absorption coefficients

- ▶ 1. Effect of Atomic Number
- ▶ 2. Effect of Density
- ▶ 3. Effect of Thickness
- ▶ 4. Effect of Gamma-Ray Energy



- ▶ the mass absorption coefficient ( $\mu_m$ ) another useful concept is defined by dividing the linear absorption coefficient by the material density

$$\mu_m = \mu / \rho \quad \dots (5)$$

Where

$\rho$  is the density of the absorber material.

## 1.4 Procedure

- ▶ 1. Connect the plugs of the electric mains.
- ▶ 2. Set the timer to 60s and the operating voltage to 380 V
- ▶ 3. Record the count rate per one minute for the background ( $I_{B.G}$ ).
- ▶ 4. Put the source in front of the GM tube.
- ▶ 5. Record the count rate ( $I_0$ ).
- ▶ 6. Place (Al) sheet midway between the source and the GM tube.
- ▶ 7. Record the count rate ( $I_1$  and  $I_2$ ) and then find  $I_{Avg}$ .

- ▶ 8. Repeat steps 6 and 7 with increasing the thickness of the absorbing material.
- ▶ 9. Place (Pb) sheet and repeat steps 6,7 and 8.
- ▶ 10. Plot a graph between  $\ln(I_0/I)$  and thickness  $X$ , if the relation is a straight line, then the absorption law is verified.
- ▶ 11. Find the slope from the graph, this is equal to the linear absorption (or Attenuation) coefficient.
- ▶ 12. Calculate the mass absorption (or Attenuation) coefficient.
- ▶ 13. Plot a graph between  $(I)$  and thickness  $(X)$ , then find the value of the half thickness graphically.
- ▶ 14. Calculate the half thickness theoretically by using eq.(4).